

# DST Series\*



GIA\* Staplers with  
Directional Staple Technology


**tyco**  

---

*Healthcare*

**Auto Suture**

# GIA\* DST Series\* Staplers



## Changing the Shape of Surgical Stapling

### GIA\* DST Staplers

REF	Description
<b>GIA6025S</b>	GIA 60-2.5 Single Use Reloadable Stapler
<b>GIA6038S</b>	GIA 60-3.8 Single Use Reloadable Stapler
<b>GIA6048S</b>	GIA 60-4.8 Single Use Reloadable Stapler
<b>SGIA6038S</b>	SGIA 60-3.8 Knifeless Single Use Reloadable Stapler
<b>GIA8038S</b>	GIA 80-3.8 Single Use Reloadable Stapler
<b>GIA8048S</b>	GIA 80-4.8 Single Use Reloadable Stapler
<b>GIA10038S</b>	GIA 100-3.8 Single Use Reloadable Stapler
<b>GIA10048S</b>	GIA 100-4.8 Single Use Reloadable Stapler

### GIA\* DST Stapler Reloads

REF	Description
<b>GIA6025L</b>	GIA 60-2.5 Single Use Loading Unit
<b>GIA6038L</b>	GIA 60-3.8 Single Use Loading Unit
<b>GIA6048L</b>	GIA 60-4.8 Single Use Loading Unit
<b>SGIA6038L</b>	SGIA 60-3.8 Knifeless Single Use Loading Unit
<b>GIA8038L</b>	GIA 80-3.8 Single Use Loading Unit
<b>GIA8048L</b>	GIA 80-4.8 Single Use Loading Unit
<b>GIA10048L</b>	GIA 100-4.8 Single Use Loading Unit
<b>GIA10038L</b>	GIA 100-3.8 Single Use Loading Unit

\*Product under development, not for sale

# Clinical Performance




Traditional Staple



**DST SERIES\*** Staple

(Staple performance comparisons based on .20 " deflection between anvil and cartridge.)

# Anvil Design




Traditional Anvil Bucket  
(Small/Narrow Target)




**DST SERIES\*** Anvil Bucket  
(Larger Target Area for Staple)

# Staple Design



Traditional Staple  
(Round Cross Section)



**DST SERIES\*** Staple  
(Square Cross Section)



## The Difference is Directional Stapling Technology



### GIA\* and TA\* SERIES\* Instruments Directional Stapling Technology

CE  
0123

[www.tycohealthcare.com](http://www.tycohealthcare.com) - [www.autosuture.com](http://www.autosuture.com)

\*Trademark ©2003. United States Surgical, a division of Tyco Healthcare Group LP. All Rights Reserved.

S-EM-GIADST/GB

**tyco**  
Healthcare

**Tyco Healthcare UK Limited**  
154 Fareham Road  
GOSPORT Hampshire  
PO 13 0AS  
UK  
Tel: +44 (0) 1329 224226  
Fax: +44 (0) 1329 224334

**Tyco Healthcare**  
Regional Export  
154 Fareham Road  
GOSPORT Hampshire  
PO 13 0AS  
UK  
Tel: +44 (0) 1329 224258  
Fax: +44 (0) 1329 224086

**Tyco Healthcare Nederland BV**  
Hogeweg 105  
NL-5301 LL ZALTBOMMEL  
The NETHERLANDS  
Tel: +31 418 57 66 00  
Fax: +31 418 57 67 91  
[TycoHealthNed@TycoHealth.com](mailto:TycoHealthNed@TycoHealth.com)  
[www.tycohealthcare.nl](http://www.tycohealthcare.nl)

**Tyco Healthcare**  
Birchwood Court Montrose Road  
P.O. Box 8108, Halfway House  
1685 MIDRAND  
SOUTH AFRICA  
Tel: +27 (0) 11 265 2000  
Fax: +27 (0) 11 265 2050

**Tyco Healthcare Belgium N.V.-S.A.**  
Generaal de Wittelaan 9/5  
B-2800 MECHELEN  
BELGIUM  
Tel: +32 (0) 15 29 44 50  
Fax: +32 (0) 15 29 44 55  
[tycohealthcare.Belgium@tyco-healthcare.com](mailto:tycohealthcare.Belgium@tyco-healthcare.com)

**Tyco Healthcare Norden AB**  
Gårdsvägen 18, Box 710  
SE-169 27 Solna  
SWEDEN  
Tel: +46-8 585 605 00  
Kundservice +46-8 585 605 05  
Fax: +46-8 585 605 06

**Tyco Healthcare Danmark**  
Langebrogade 6E, 4.th.  
DK-1411 København K  
DENMARK  
Tel: +45-70 27 52 50  
Kundeservice +45-70 27 52 50  
Fax: +45-70 27 51 50

**Tyco Healthcare Norge**  
Løkketangen 20  
N-1337 Sandvika  
NORWAY  
Tel: +47-67 80 55 10  
Kundeservice +47-67 80 55 18  
Fax: +47-67 80 55 11

**Tyco Healthcare Finland**  
Köydenpunojankatu 2D, 6 krs.  
PL 340  
FIN-00181 Helsinki  
FINLAND  
Puh +358-(0)9-6226 840  
Asiakaspalvelu  
+358-(0)9-6226 8418  
Fax: +358-(0)9-6226 8460